

Medical Portal System – API Documentation

1. Patient Login and View Reports

POST /api/auth/login

Authenticate a patient and return a session token.

Request:

```
{
  "email": "patient@example.com",
  "password": "mypassword"
}
```

Response (200):

```
{
  "status": "success",
  "token": "jwt_token_here",
  "userId": "patient_123"
}
```

Errors:

- 401 → Invalid credentials
- 403 → Account locked

GET /api/patients/{patientId}/reports

Fetch all reports for a patient.

Headers: Authorization: Bearer <token>

Response (200):

```
{
  "reports": [
    {
      "reportId": "rpt_001",
      "date": "2025-10-01",
      "riskScore": 0.82,
      "riskLevel": "HIGH"
    }
  ]
}
```

```
]
}
```

GET /api/reports/{reportId}
Get detailed report data including AI analysis.

Headers: Authorization: Bearer <token>

Response (200):

```
{
  "reportId": "rpt_001",
  "patientId": "patient_123",
  "imageUrl": "/storage/xray_123.png",
  "aiAnalysis": {
    "riskScore": 0.82,
    "confidence": 0.91,
    "recommendations": ["Schedule CT scan", "Consult specialist"]
  }
}
```

2. Radiologist Upload X-Ray Report

POST /api/reports/upload
Upload a patient X-Ray for AI analysis.

Form-Data Request:
patientId=patient_123
file=xray_image.png

Response (200):

```
{
  "reportId": "rpt_456",
  "status": "processing",
  "riskScore": null
}
```

Errors:

- 400 → Invalid image
 - 413 → File too large
-

POST /api/ai/analyze

Run AI analysis on the uploaded image.

Request:

```
{
  "reportId": "rpt_456",
  "imageUrl": "/storage/tmp/xray_123.png"
}
```

Response (200):

```
{
  "riskScore": 0.65,
  "riskLevel": "MEDIUM",
  "confidence": 0.88,
  "analysisSummary": "Possible pneumonia patterns detected."
}
```

POST /api/notifications/send

Notify stakeholders about report and risk results.

Request:

```
{
  "reportId": "rpt_456",
  "recipients": ["patient_123", "doctor_001"],
  "type": "RISK_ALERT",
  "riskLevel": "MEDIUM"
}
```

Response (200):

```
{ "status": "sent" }
```

3. Doctor Review AI-Generated Report

GET /api/reports/priority?doctorId={doctorId}

Get prioritized report queue for a doctor.

Response (200):

```
{
  "reports": [
```

```
{ "reportId": "rpt_456", "riskLevel": "HIGH", "status": "pending" }
]
```

GET /api/reports/{reportId}/details
Get report details and AI analysis.

Response (200):

```
{
  "reportId": "rpt_456",
  "imageUrl": "/storage/xray.png",
  "aiAnalysis": {
    "riskScore": 0.65,
    "recommendations": ["Follow-up imaging"]
  }
}
```

POST /api/reports/{reportId}/review
Submit a doctor's review of an AI-generated report.

Request:

```
{
  "doctorId": "doctor_001",
  "notes": "Confirm pneumonia, suggest antibiotics.",
  "urgencyLevel": "HIGH",
  "finalDiagnosis": "Pneumonia"
}
```

Response (200):

```
{ "status": "reviewed", "reviewId": "rev_001" }
```

4. Tech Team Monitor AI Model Performance

GET /api/models/active
Fetch currently deployed models.

Response (200):

```
{
  "models": [
```

```
{ "modelId": "v3.2", "status": "active", "deployedAt": "2025-09-15" }
]
```

GET /api/models/{modelId}/metrics
Fetch performance metrics for a model.

Response (200):

```
{
  "accuracy": 0.92,
  "precision": 0.88,
  "recall": 0.90,
  "f1Score": 0.89,
  "timeRange": "last_7_days"
}
```

POST /api/models/{modelId}/retrain
Trigger model retraining.

Request:

```
{ "dataset": "lung_xray_2025", "epochs": 100 }
```

Response (200):

```
{ "status": "training_started", "jobId": "train_789" }
```

POST /api/models/{modelId}/deploy
Deploy a new trained model.

Request:

```
{ "version": "v3.3", "jobId": "train_789" }
```

Response (200):

```
{ "status": "deployed", "modelId": "v3.3" }
```

5. Complete User Journey (End-to-End)

Key APIs used:

1. Radiologist → POST /api/reports/upload

2. System → POST /api/ai/analyze
3. Notification → POST /api/notifications/send
4. Patient → POST /api/auth/login → GET /api/patients/{id}/reports
5. Doctor → GET /api/reports/priority → POST /api/reports/{id}/review
6. Patient follow-up → POST /api/appointments/schedule

6. AI Model Prediction Pipeline

POST /api/ai/validateImage
Validate uploaded image.

Response:
{ "valid": true, "issues": [] }

POST /api/ai/preprocess
Preprocess image before AI analysis.

Response:
{ "processedImageUrl": "/storage/tmp/processed_xray.png" }

POST /api/ai/extractFeatures
Extract features for ML model.

Response:
{ "featureVector": [0.12, 0.55, 0.87, ...] }

POST /api/ai/predict
Run forward pass prediction.

Response:
{ "rawPredictions": [0.1, 0.3, 0.6] }

POST /api/ai/calculateRisk
Calculate risk score from predictions.

Response:
{

```
"riskScore": 0.6,  
"riskLevel": "MEDIUM",  
"explanation": "Detected anomalies in left lung."  
}
```

POST /api/reports/{reportId}/saveAnalysis
Save AI analysis results.

Response:

```
{ "status": "saved", "analysisId": "anl_123" }
```